



DEVICE DETAILS

NAME OF DEVICE	KF99-SURGICAL FACE MASK
ESTABLISHMENT NAME	ICON PACKAGING SDN BHD
REGISTRATION NO	GA7339622-91110
BRAND NAME	CHAN.GER
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A single use non sterile and latex free device made of non-woven, structural filter and melt blown. Adjustable to fit faces of various shapes and sizes. The face mask tested for Bacterial Filtration Efficiency (BFE), Particle Filtration Efficiency (PFE) and Viral Filtration efficiency (VFE) that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight.
DEVICE INTENDED PURPOSE	A mask is meant to help block large particle droplets, splashes or splatter that may contain germs (viruses & bacterial) keeping it from it reaching to your mouth. Face mask may also help to reduce the exposure of your saliva and respiratory secretion to others. This device, is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
VALIDITY DATE OF REGISTRATION	25/04/2022 - 24/04/2027

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	KF99 SURGICAL FACE MASK	CG7001
2	KF99 SURGICAL FACE MASK	CG7006