



DEVICE DETAILS

NAME OF DEVICE	MOVEN Lightweight Transport Wheelchair
ESTABLISHMENT NAME	BEST MADE PHARMACEUTICAL SDN BHD
REGISTRATION NO	GA9618422-91832
BRAND NAME	Lightweight Transport Wheelchair / MOVEN
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Lightweight wheelchair made of aluminum frame, foldable backrest, fixed armrest, detachable legrest, double crossbar and 8"/12" front/back wheels.
DEVICE INTENDED PURPOSE	To help individuals who have impairments that limit their ability to walk to move.
VALIDITY DATE OF REGISTRATION	07/05/2022 - 06/05/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	MOVEN Lightweight Transport Wheelchair	BM01LABJ-12



MEDICAL DEVICE AUTHORITY