



DEVICE DETAILS

NAME OF DEVICE	Surgical Face Mask
ESTABLISHMENT NAME	U CARE PROTECTION INDUSTRIES SDN BHD
REGISTRATION NO	GA3161022-94988
BRAND NAME	Surgical Face Mask / Daxxo
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	A single use non sterile and latex free device made of 2 layer of non woven and 1 layer of melt blown. Adjustable to fit faces of various shapes and sizes. The middle layer of face mask (Melt blown) which tested for Bacterial Filtration Efficiency that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight. This device is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
DEVICE INTENDED PURPOSE	A disposable single use surgical face mask is to be used as a barrier to prevent air-borne transfer of infectious, large particle droplets and various micro-organisms. It is designed to be used in hospital, clinic, operating theatre, dental, cleanroom and critical environment.
VALIDITY DATE OF REGISTRATION	20/06/2022 - 19/06/2027

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	Surgical Face Mask	DX4P101