



DEVICE DETAILS

NAME OF DEVICE	Medical Face Mask 3 ply
ESTABLISHMENT NAME	THE BEST PERFECT SDN. BHD.
REGISTRATION NO	GA5556722-83996
BRAND NAME	Nine-Care
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A single use non sterile and latex free device made of 2 layer of non woven and 1 layer of melt blown. Adjustable to fit faces of various shapes and sizes. The middle layer of face mask (Melt blown) which tested for Bacterial Filtration Efficiency that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight. This device is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
DEVICE INTENDED PURPOSE	A surgical face mask is meant to help block large particle droplets, splashes or splatter that may contain germs (viruses & bacteria), keeping it from reaching to your mouth. Face mask may also help reduce exposure of your saliva and respiratory secretions to others.
VALIDITY DATE OF REGISTRATION	24/01/2022 - 23/01/2027

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	3PLY DISPOSABLE MEDICAL FACE MASK	Earloop
2	3PLY DISPOSABLE MEDICAL FACE MASK	HIJAB Headloop