



DEVICE DETAILS

NAME OF DEVICE	KN95 MEDICAL FACE MASK
ESTABLISHMENT NAME	TNC RETAILS ENTERPRISE
REGISTRATION NO	GA5201222-106084
BRAND NAME	CAREION X MASKING
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A single use non sterile and latex free device made of non-woven, structural filter and melt blown. Adjustable to fit faces of various shapes and sizes. The face mask tested for Bacterial Filtration Efficiency (BFE) that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight.
DEVICE INTENDED PURPOSE	A mask is meant to help block large particle droplets, splashes or splatter that may contain germs (bacterial) keeping it from it reaching to your mouth. Face mask may also help to reduce the exposure of your saliva and respiratory secretion to others. This device, is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
VALIDITY DATE OF REGISTRATION	14/10/2022 - 13/10/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KN95 MEDICAL FACE MASK	KN95/ADULT
2	KN95 MEDICAL FACE MASK	KN95/KID