



DEVICE DETAILS

NAME OF DEVICE	FINGERTIP PULSE OXIMETER
ESTABLISHMENT NAME	KINSMEDIC SDN BHD
REGISTRATION NO	GB2362522-109761
BRAND NAME	CHOICEMMED
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO2. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO2 level.
DEVICE INTENDED PURPOSE	Fingertip pulse oximeter is a handheld non-invasive device intended for spot checking of oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult and pediatric patients in hospitals, hospital type facilities and homecare.
VALIDITY DATE OF REGISTRATION	22/11/2022 - 21/11/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		