



DEVICE DETAILS

NAME OF DEVICE	Non-woven Face Mask
ESTABLISHMENT NAME	NUTRIMEDIQ SDN BHD
REGISTRATION NO	GA10783520-38976
BRAND NAME	Dermaseal Non-woven Face Mask
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	3 ply is made of Polypropylene Spunbond and Meltblown nonwoven fabrics with adjustable nose strip and flat or round elastic ear- loop to give basic protection to the user.
DEVICE INTENDED PURPOSE	Basic protection from airborne particle and bacteria. Mainly used in medical, dental, cleanroom, critical environment, food handling/production and industrial manufacturing.
VALIDITY DATE OF REGISTRATION	17/01/2020 - 16/01/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Non-woven Face Mask	9-555391-821518

