



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER OLV-4201
ESTABLISHMENT NAME	NIHON KOHDEN MALAYSIA SDN BHD
REGISTRATION NO	GC95939699418
BRAND NAME	NIHON KOHDEN (OLV-4201)
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	THE PULSE OXIMETER IS A COMPACT MONITOR WITH A 7-INCH DISPLAY AND IS DESIGNED FOR TOUCHSCREEN OPERATION AT THE PATIENT'S SIDE. THE PULSE OXIMETER IS INSTALLED NEAR THE PATIENT TO DISPLAY THE PATIENT'S VITAL SIGNS (SPO2, PULSE RATE, ETC.) ON THE DISPLAY AND GENERATE ALARMS. NEONATES, CHILDREN AND ADULTS CAN BE MONITORED WITH THIS PULSE OXIMETER. WHEN THE OPERATION MODE IS SET TO NETWORK MODE, THE PULSE OXIMETER CAN CONNECT TO A NIHON KOHDEN CENTRAL MONITOR NETWORK AND COMMUNICATE WITH CENTRAL MONITORS AND BEDSIDE MONITORS IN THE NETWORK.
DEVICE INTENDED PURPOSE	THE PULSE OXIMETER IS INSTALLED NEAR THE PATIENT TO DISPLAY THE PATIENT'S VITAL SIGNS (SPO2, PULSE RATE, ETC.) ON THE DISPLAY AND GENERATE ALARMS. THIS PULSE OXIMETER MUST BE USED BY QUALIFIED MEDICAL PERSONNEL WITH FULL KNOWLEDGE OF ITS FUNCTIONS AND OPERATIONS WITHIN THE ENVIRONMENT OF HOSPITALS SUCH AS CRITICAL CARE, ICU, HCU, OR, GCU, OUTPATIENT CLINIC, GENERAL WARD AND SO ON. IT IS NOT INTENDED FOR PATIENTS TO USE. IT IS NOT INTENDED TO USE OUTSIDE OR AT HOME, EITHER.
VALIDITY DATE OF REGISTRATION	22/02/2023 - 21/02/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		