



DEVICE DETAILS

| | |
|-------------------------------|---|
| NAME OF DEVICE | FRONTIFIE 4-PLY SURGICAL FACE MASK |
| ESTABLISHMENT NAME | NEUTROVIS SDN BHD |
| REGISTRATION NO | GA3267123-132590 |
| BRAND NAME | FRONTIFIE |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | This single use surgical face mask with BFE \geq 98% and PFE \geq 98% is made from premium medical grade melt-blown filter fabric that will protect you from germs, bacteria and viruses. The soft, non-woven, inner material of the mask provides better comfort and reduced irritation. It is able to filter out live bacteria (viable particles) from 5 to 1 micron and micro particles (nonviable particles) from 1 to 0.1 micron |
| DEVICE INTENDED PURPOSE | 4 PLY DISPOSABLE MEDICAL AND SURGICAL FACE MASK IS A NON-WOVEN MATERIAL THAT OFFERS FILTRATION-EFFICIENT, FLUID-RESISTANT, SUPERIOR COMFORT AND BREATHABILITY TO HELP REDUCE WEARER TO EXPOSURE OF POTENTIALLY INFECTIOUS PARTICLES, BLOOD AND BODY FLUIDS . IT IS DESIGNED TO BE USED IN HOSPITAL ESPECIALLY IN OPERATING THEATER AND CLEAN ROOM ENVIRONMENT. |
| VALIDITY DATE OF REGISTRATION | 26/05/2023 - 25/05/2028 |

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

| NO | NAME OF DEVICE | IDENTIFIER |
|----|---|------------|
| 1 | Frontifie 4-Ply Surgical Face Mask - Admiral Blue | FTF003 |
| 2 | Frontifie 4-Ply Surgical Face Mask - Bold Black | FTF002 |