



DEVICE DETAILS

NAME OF DEVICE	HERA W9 DIAGNOSTIC ULTRASOUND APPARATUS
ESTABLISHMENT NAME	LAC MEDICAL SUPPLIES SDN BHD
REGISTRATION NO	GB2791320-44450
BRAND NAME	HERA W9 - SAMSUNG MEDISON
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	<p>The HERA W9 is general purpose, mobile, software controlled, diagnostic ultrasound system. Its function is to acquire ultrasound data and to display the data as B-mode, M-mode, Pulsed wave (PW) Doppler, Continuous wave (CW) Doppler, Color Doppler, Tissue Doppler Imaging (TDI), Tissue Doppler Wave (TDW), Power Amplitude Doppler, Pulse Inversion Harmonic Imaging (S-Harmonic), Directional Power Doppler (S-Flow), Color M-Mode, 3D Imaging Mode, 4D Imaging Mode, Elastoscans+ Mode, Tissue Harmonic Imaging, MV-Flow Mode or as a combination of these modes. The HERA W10 also gives the operator the ability to measure anatomical structures and offers analysis packages that provide information that is used to make a diagnosis by competent health care professionals. The HERA W9 has real time acoustic output display with two basic indices, a mechanical index and a thermal index, which are both automatically displayed.</p>
DEVICE INTENDED PURPOSE	<p>The ultrasound diagnostic system and probes are designed to obtain ultrasound images and analyze body fluid. The clinical applications include: Fetal/Obstetrics, Abdominal, Gynecology, Pediatric, Small Organ, Neonatal Cephalic, Adult Cephalic, Trans-rectal, Trans-vaginal, Muscular-Skeletal (Conventional, Superficial), Urology, Cardiac Adult, Cardiac Pediatric and Peripheral vessel.</p>
VALIDITY DATE OF REGISTRATION	24/06/2020 - 23/06/2025

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		